

A Fast and Sensitive Method for Residual Hydrazine Analysis in Pharmaceutical Samples

Vera Leshchinskaya, Priscilla Richberg, Yury Zelechonok, and Anne Kelly
Bristol Myers Squibb Company, P.O.BOX 4000, Princeton, NJ, 08543,*

** SIELC Technologies, 65 E. Palatine Rd., Suite 221, Prospect Heights, IL 60070*

Hydrazine Properties



Hydrazine is

- a highly reactive molecule widely used in chemical syntheses of intermediates and active pharmaceutical ingredients
- a known carcinogen - can have an impact on product risk assessment if present in the final drug substance and drug product as an impurity.
- a small, polar, basic molecule with no chromophore



Objective

Develop a method for detection of low (ppm) level of hydrazine in pharmaceutical samples and wash solutions. Method criteria:

- simple sample preparation step
- efficiency
- sensitivity
- ease of transfer



Control of Genotoxic Impurities in Pharmaceuticals

According to current regulatory guidance, it is assumed that genotoxic compounds have the potential to damage DNA at very low level of exposure. Thus, actual and potential impurities most likely to arise during synthesis, purification, and storage should be identified. Both the EMEA guidelines and a PhRMA white paper propose a limit of 1.5 µg/day for genotoxic impurities in pharmaceuticals.



Analytical Challenges in Setting Limits for Genotoxic Impurities

- Low level (ppm) needs to be detected
- Many genotoxic impurities are highly reactive and therefore difficult to analyze
- Sample matrix can be chemically similar and present at extremely high level
- Techniques being employed aren't standard and may present significant technology transfer issues




Analytical Background

Existing methods for hydrazine

- are mostly applicable to environmental samples (water, air, soil)
- employ sophisticated equipment and/or derivatization
- tend to be cumbersome and time consuming
- are subject to interference, especially by hydrazine derivatives (colorimetric and titrimetric methods)



GC Analysis of Residual Hydrazine




Advantages of using Acetone:

- Serves as both diluent and derivatizing agent
- Symmetrical molecule (produces only one isomer)
- Fast reaction rate
- Low toxicity
- Good GC separation from the product of the reaction.




GC Analysis of Acetone Azine



Column: 20-m length x 0.18-mm i.d., 1 μ m film thickness, 6% cyanopropylphenyl-94% dimethylpolysiloxane (DB-624);
Detector: FID, 280°C; Injector: 1 μ L into the split port maintained at 200 °C;
Carrier Gas: Helium; Column Flow: 1.3 ml/min; Split ratio/Split Flow: 110/143 ml/min.

Linearity Data for Acetone Azine





Concentration		Peak Area
mg/mL	ppm	
0.000842	0.84	0.1168
0.005052	5.05	0.5077
0.00842	8.40	0.8384
0.05052	50.5	5.2658
0.0842	84.2	8.8097
0.5052	505.2	54.1033
0.842	842.0	88.3914


S/N at 0.84 ppm is ~10:1
0.84 ppm of Acetone azine
corresponds to 0.24 ppm
of free hydrazine in solution
or 2.4 ppm/100 mg of sample



Conversion of Hydrazine Mono-Hydrate to Acetone Azine: GC/MS Data



Conversion Rate of Hydrazine Mono-Hydrate to Acetone Azine



Time, min	% Conversion
5	10.3
103	22.8
153	29.1
170	31.7
187	34.1
204	36.8
259	37.5
276	39.7
573	69.1
590	77.0



Reaction Conversion of Hydrazine Mono-Hydrate to Acetone Azine Catalyzed by Acid

Time, min	% Conversion, Acetone Hydrazone/ Acetone Azine	
	HCOOH	CH ₃ COOH
0	0/100	5.5/94.5
15	0/100	0/100

Conclusion : The derivatization proceeds in less than 3 min in the presence of formic acid and within ~15 min with Acetic Acid.



Reaction Conversion of Hydrazine Dihydrochloride to Acetone Azine

■ Problem:

Low solubility of the salt in acetone

Conversion rate is ~50% within 24 hours

■ Solution:

Standard prep

- dissolved~ 10 mg/2 ml of H₂O
- diluted with acetone to 200 ml

The derivatization proceeds in 2.5 hours (97.5%). Presence of formic acid did not show improvement in the reaction rate



Sample Analysis for Residual Hydrazine.

■ Problem:

~100 mg needed to be dissolved in 1 ml of acetone for ppm level hydrazine analysis.

■ Solution:


Sample prep

- added 1 ml of 1/99 H₂O/Acetone to 100 mg of the sample.
- Heated at ~ 40°C. Mixed well.

The solution became clear for a short period of time (1-2 minutes) and then the precipitate formed again under cooling to room temperature.

■ Question:


Did the reaction occur or did the recrystallization take place?



Compound A




^1H NMR Spectra Comparison (1)



DMSO- d_6 , 400 MHz

Compound A



Precipitate from Acetone,
yield is $\sim 100\%$




Bristol-Myers Squibb Company


^1H NMR Spectra Comparison (2)

DMSO-d₆, 400 MHz

Acetone layer



Precipitate from Acetone spiked
with Acetone Azine (AA)



Bristol-Myers Squibb Company

NMR Study of reaction conversion

- NMR of the precipitate showed mostly hydrazone, the conversion is close to 100% by weight
- NMR of supernatant showed ~30% of unreacted sample, hydrazone and acetone azine formed due to the presence of residual hydrazine in the sample
- Spiking experiment showed 100.9% recovery of spiked hydrazine (synergistic effect?)



Residual Hydrazine Analysis in Neutral Molecules with Low Solubility in Acetone



Procedure:

- Dissolved sample (~50mg) in NMP or any other suitable solvent
- Added 0.1 ml of acetone/HCOOH;
- Analyzed sample by the method described above

The recovery of hydrazine from solution fortified with hydrazine mono-hydrate was 103%

Direct Analysis of Hydrazine in Aqueous Solutions

Mixed-Mode Chromatography:




Column: Primesep 100 (Sielc Technology)

Gradient: 20%ACN/80%H₂O/0.2%TFA for 6 min, then 80%ACN/20%H₂O/0.2%TFA, 1 ml/min, Sample concentration - 31.4 mg/ml, inj. Vol.- 10 μ l

Detection: ELSD; T= 40°C; Photomultiplier voltage- 500


Direct Analysis of Hydrazine by Mixed-Mode Chromatography

ELSD



Detection limit - ~100 ppm

CAD



Detection limit - ~35 ppm

Conclusion

- The reaction of acetone with hydrazine free base and HCl salts has been investigated and conditions were developed for quantitative and fast conversion to acetone azine
- A fast and sensitive GC method for residual hydrazine analysis in pharmaceutical samples has been developed
- Excellent linear relationship between acetone azine concentration and FID response was demonstrated with minimal detectable concentration 0.24 ppm
- The method is applicable to pharmaceutical compounds at all stages of development with some modifications
- Hydrazine can be analyzed in aqueous solutions using mixed-mode chromatography and CAD detection with minimal detectable concentration 35 ppm

Acknowledgments

- Su Pan (CAD data)
- Kana Yamamoto
- Reginald Cann
- Xinhua Qian