


SIELC Technologies, Inc.
Wheeling, IL 60090 USA
P. 847-229-2629 F. 847-655-6079
mail@sielc.com www.sielc.com

SIELC

MEASURING EDTA USING HPLC



EDTA Analysis Background

Ethylenediaminetetraacetic acid (EDTA) has multiple industrial, medical, biological and research applications. It presents a number of difficult problems when one attempts to analyse it using HPLC:

- It has no retention in reversed-phase chromatography
- It possesses multiple charges at acidic and neutral pH, which are most common in HPLC
- EDTA molecules have low UV absorption above 230 nm
- It has strong chelating properties and gets retained by all metal components of the HPLC system
- It is usually a part of the formulation, mixed with surfactants and other components, which can complicate separation and detection

SIELC developed a new selective methodology to address these difficulties and offers a simple and reliable method for EDTA quantitation in most liquid samples. The method is based on forming a complex of EDTA with Fe(III), which has strong UV absorption maximums at 260 nm (Fig. 1).

This complex can be measured by a UV-Vis detector and can be separated from the signal of complex-free Fe(III) using Newcrom BH HPLC column (Fig. 2).




Fig. 1. UV spectra of EDTA • Fe(III) complex, FeCl₃ solution, and EDTA solution.

EDTA Analysis

EDTA Standards Solution A

For the preparation of EDTA standard solution, 5 mg of EDTA were accurately weighed, transferred into a 5 mL volumetric flask and dissolved in 0.001M NaOH aqueous solution using sonication or magnetic stirrer. The EDTA stock solution (1.0 mg/mL) should be stored in a cold dark place and can be used for a week to prepare standards of required concentration.




Fig. 2. Chromatogram of EDTA complex with Fe(III). HPLC column: Newcrom BH, 4.6 x 150 mm, 5 μ m. Flow rate: 1.0 mL/min. Mobile phase: 2% MeCN in water containing 0.1% H₂SO₄. Detection: UV 260 nm.




Fig. 3. Isocratic chromatogram of EDTA complex with Fe(III). HPLC column: Newcrom BH, 4.6 x 150 mm, 5 μ m. Flow rate: 1.0 mL/min. Mobile phase: 2% MeCN in water containing 0.1% H₂SO₄. Detection: UV 260 nm.

Iron(III) chloride Solution B

The standard stock solution of Iron(III) chloride (10 mg/mL) was prepared in water. 50 mg of Iron(III) chloride hexahydrate were accurately weighed, transferred into a 5 mL volumetric flask, dissolved in water and sonicated, if needed.

General procedure for Iron(III) - EDTA complex analysis

Mix 100 μ L of Solution A (or unknown sample) with 300 μ L of Solution B and 600 μ L of water. Place this mixture in a plastic HPLC vial for analysis. HPLC conditions described in Fig. 2.

For situations when a sample has more than one component that needs to be analyzed, gradient conditions can be used as well, in the same run. Example below (Fig. 3) shows isocratic conditions for EDTA testing only. The next example (Fig. 4) shows gradient run when other ingredients of the formulation can be measured at the same time.





Fig. 4. Gradient method. HPLC column: Newcrom BH, 4.6 x 150 mm, 5 μ m. Flow rate: 1.0 mL/min. Mobile phase: MeCN gradient from 2 to 70% in water in 10 min, mobile phase contains 0.1% H₂SO₄. Detection: UV 260 nm. Injection: 5 μ L of sample.

Summary

Iron(III)-EDTA complex response vs concentration
Range 1 - 10 µg/mL



Iron(III)-EDTA complex response vs concentration
Range 10 - 100 µg/mL




Fig. 5. Linearity study of the EDTA analysis quantitation method. Analytical column: Newcrom BH, 4.6 x 150 mm, 5 µm. Flow rate: 1.0 mL/min. Mobile phase: 2% MeCN in water containing 0.1% H₂SO₄. Detection: UV 260 nm. Injection: 5 µL of EDTA standards with FeCl₃ solution.

SIELC




Fig. 6. Schematic structure of Newcrom BH column surface chemistry

Linearity of this method was investigated using increased concentrations of EDTA. Good linearity of the UV absorption response at 260 nm of EDTA complex with Fe(III) ions was observed within a wide concentration range.

Simple sample preparation includes mixing the unknown sample with Fe(III) stock solution followed by an HPLC separation. Sensitivity (LOQ) of down to 0.1 ppm of EDTA in samples has been achieved routinely.

Column surface chemistry is schematically presented in Fig. 6. It is comprised of a long alkyl chain with a terminal strong basic group carrying a permanent positive charge. This ligand is covalently bonded to the silica gel support. Silica gel employed is of spherical type, with fully porous particles of 5 µm in diameter and pore size of 120 Å.



Column part number
NBH-150.46.0510

To order a column or ask a question
send your message to sales@sielc.com
or call us at +1 (847) 229-2629

SIELC Technologies, Inc.

www.sielc.com email: mail@sielc.com ph. 847-229-2629 fax 847-655-6079