



# Betula Species HPLC Comparison.

**Column:** Primesep B2  
**Size:** 150 x 4.6 mm  
**Mobile phase:** MeCN gradient  
from 70 to 90 in 15 min  
with 0.05% H<sub>2</sub>SO<sub>4</sub>  
**Flow:** 1.0 ml/min  
**Detection:** UV 210 nm

**Column:** Primesep 100  
**Size:** 150 x 4.6 mm  
**Mobile phase:** MeCN gradient  
from 70 to 90 in 15 min  
**Flow:** 1.0 ml/min  
**Detection:** UV 210 nm

